



**Proudly Made in Malaysia**  
Protecting Hands Since 1988



**Careplus**

Quality and Relationship You Can Trust



**For Every Precise Movement**

## Nitrile Gloves

### Examination Grade Powder Free

A Nitrile glove provides protection from common chemicals, with good puncture resistance and tear strength



Single-use



Latex free



Powder free

## Product Specification

Meet or exceed **ASTM D6319 & EN 455**

- Color: Blue, White, Indigo Blue, Violet, Black
- Weight: Available in 3.5gm, 4.6gm and 6gm (standard cuff)
- Length: 240mm
- Sizes: XS, S, M, L, XL
- Packing: 100 gloves per dispenser, 10 dispensers per carton
- \*Other packing options are available*



Blue color



White color



Black color



Violet color



Indigo Blue color



Dark Blue color

## Features

- **Excellent fit and soft feel** for outstanding comfort and protection
- **Latex Free** for users who suffer from natural rubber latex allergies
- **Textured finger tip** enhances grip in wet and dry conditions
- **Nitrile provides protection** from common chemicals, with good puncture resistance and tear strength



All pictures shown are for illustration purpose only



**Careplus (M) Sdn Bhd**

(Fully owned subsidiary of the public listed company, Careplus Group Berhad)

Lot 120 & 121, Jalan Senawang 3, Senawang Industrial Estate, 70450 Seremban, Negeri Sembilan Darul Khusus, MALAYSIA.



info@careplus.com



+60 6677 2781



www.careplus.com



Careplus Malaysia